
de	Gebrauchsanweisung: Stimmgabeln.....	2
en	User's Manual: Tuning Forks	5
fr	Mode d'emploi : Diapasons	8
it	Istruzioni per l'uso: Diapason	11
es	Instrucciones de empleo: Diapasones.....	14
pt	Manual de operação: Diapasões	17
ru	Руководство по применению: Камертоны ..	20

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
User's Manual

Rydel Seiffer, Lucae and Hartmann Tuning Forks

Example of Rydel Seiffer tuning fork,
c 128/C 64 Hz with damper and
plastic foot

adjustable

Example of Lucae tuning fork,
c2 512 Hz without dampers
with plastic foot

Example of Hartmann tuning fork,
c2 512 Hz without damper
and plastic foot

Dear buyer! Thank you for choosing a KaWe product. Our products are known for high quality and durability. This KaWe device complies with the provisions of Regulation EU 2017/745 (European Regulation on Medical Devices), whereby it is classified as a medical product of class 1.

 **Please carefully consider this instruction and read all the provisions prior to use. Follow the maintenance recommendations.**

Intended use: KaWe tuning forks are used as a source of sound conduction through the air and through the bones. They are provided for medical purposes, for example, hearing tests, sensitivity tests and other diagnostic purposes.

Group of patients with indications for use: All

Application: Only qualified medical professionals can use the tuning fork. The tuning fork is used as a source of sound conduction through the air and through the bones. Use on healthy skin only. The vibrating fork is placed on a plastic stand or on the studied areas of the body with its end. First of all, these areas include places on the limbs where the joints are not hidden by muscles or adipose tissue. To ensure long-term control, compare the data at the same points of the study. This requires the patient telling about the point of the fork attenuation when vibration is still felt.

Application according to Rydel-Seiffer: The Rydel-Seiffer tuning fork produces two different frequencies: with the silencers installed, the tone C = 64 Hz (as noted on the silencers) and the tone c = 128 Hz with the silencers removed (as indicated on the tuning fork itself).

If the silencers are on, their lower edge should be installed on the mark of the tuning fork teeth. At the same time, the marking "c 64" on the silencers and the marking "c 128" on the tuning fork (front),

as well as the label, must always be turned to the observer. The silencers can be fixed and removed using thumb screws. Tighten the thumb screws firmly.

Scales are applied on the silencers on the reverse side of the tuning fork. Read the intensity of the tuning fork vibrations when examining the patient with the same vibration force by observing these scales. This happens as follows:

- A tuning fork with installed silencers (i.e., with a tone of C = 64 Hz) hits the ball of a thumb. The tuning fork teeth vibrate strongly at first, and both triangles are completely blurred. After a few seconds, depending on the force of the impact, the main line triangles become clearer again and increase slowly, i.e. with a slowly fading intensity of vibrations to the first, second, third, etc. lines on the scale. Thus, the examination can begin with a certain mark on the scale. Depending on the incidence of light or better visibility, a black or white triangle can be observed.

It is possible to draw a conclusion in advance by weakened vibration sensation about the initial damage in the nervous system. A reproducible vibration value is a precondition for the practical application of this knowledge.

Application as per Lucae: The Lucae tuning fork produces 7 different frequencies, from c 128 to h 240. The exact frequency values can be found in the following table.

	c	d	e	f	g	a	h
Frequency, Hz	128	144	160	170,66	192	213,33	240

The frequency changes by shifting the weights on the teeth. These weights can be rearranged manually. When changing the pitch, make sure that the spring in the weight enters the groove of the tuning fork. This is the only way to achieve an accurate frequency. Both weights should be installed at the same height. Therefore, medium pitch cannot be generated by setting one weight to one pitch, and the second to another. The stamped letter right under the weight indicates the generated pitch.

Tuning fork should hit on the ball of the thumb, on wood or plastic, but in no case on solid objects (for example, an iron sheet, steel or stone). The beat should be short and elastic, i.e. the touch with the tuning fork should be very short. The pitch oscillation would be the longest due to such a short beat.

During the examination (sensitivity test, etc.), the tuning fork is placed after the beat with the base of the black stand on the examined points of the patient.

Improper use of the device/contraindications for use: Any other use of use outside the above intended use shall be considered improper. The manufacturer is not responsible for the damage caused as a result of such use; the responsibility for the risks in this case rests solely with the user.


Additional information, maintenance and storage: The device will serve for years with proper handling and storage.

⚠ Warnings: Please note that minor damage may affect the frequency of vibrations. If the tuning fork is equipped with screwed weights, they can not be unscrewed. The frequency is changed by moving the weights. Only the manufacturer can perform the exact adjustment.

Secondary treatment limitations: Secondary treatment has a minor effect on these devices. The service life of the product is usually defined based on the degree of wear and use-related damage. Prompt secondary processing of the tool is recommended after its application.

Instructions

Transportation and storage: Rydel-Seiffer: Tuning the fork after transportation: The silencers on the tuning forks may get loose during transportation. Before using, make sure that both silencers are fixed to the tuning fork so that the marking on the back of the fork is located under the lower edge of the silencers (see figure). Only in this case the correct application of the fork (frequency) may be guaranteed.

Preparation for cleaning: No special requirements.

Manual cleaning: Do not clean the device in sterilizers after using. Wipe the tuning fork and the stand with a dry or slightly oily cloth to clean them.

Automatic cleaning: Make sure to use a cleaning machine that is designed for medical devices. Use only cleaning products suitable for metals. Observe the concentrations and exposure time specified by the manufacturer.

Disinfection: The parts in contact with the patient should be disinfected after each use. Use alcohol-based disinfectants; observe the data and exposure time specified by the manufacturer.

Maintenance: Remove or replace the damaged parts.

Control and functional check: Check the ease of movement of moving parts, for example, silencers. Visually check for damages and wear. Before using, check the thread tightness.

Package: Separate. A standard plastic bag can be used. The bag must be of sufficient size for the tool so that the clasp is not under tension.

Sterilization: It is a non-surgical device; hence, sterilization is not required or possible.

Storage: Keep the tool away from dust, humidity and contamination during storage and transportation.

Ambient temperature: Room temperature:

Disposal: The product can be disposed of to the scrap metal collection point.

General accessories: To find out more about the products, visit our homepage: www.kawemed.com.

Manufacturer: **KaWe**

Manufacturer's contact details: Dealer's address or phone number: +49-7141-68188-0.

Basic UDI-DI: 4030155KaWe08SY

Warranty: The warranty for the device totals two years from the date of purchase with proper handling and following the operation manual. If you have any additional questions or need repair, please, contact your dealer.

Symbol key:

	Attention
	Manufacturer & date of manufacture
	Lot code
	Please follow the User's Manual
	CE conformity mark
	Article number
	Medical device
	UDI Data Carrier

Information for users and patients

All serious incidents involving this device must be immediately reported to the manufacturer and the competent authorities of the respective Member State in which the user and/or patient is located.

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KIRCHNER & WILHELM GmbH + Co. KG