GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403
Export: tel. +39 02 953854209/221/225 - fax +39 02 95380056
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

SCHIZZETTONI EAR WASHING SYRINGES SERINGUES POUR LAVAGES AURICULAIRES SPRITZEN ROCIADOR IRRIGADORES ΨΕΚΑΣΤΗΡΕ أدوات لرش الأدوية

Manuale d'uso - User manual
Manuel de l'utilisateur
Gebrauchs- und Instandhaltungsanleitung
Guia de uso - Guia para utilização
Οδηγίες χρήσης - دليل الإستعمال والرعاية

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.
ATTENTION: The operators must carefully read and completely understand the present manual before using the product.
AVIS: Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.
ACHTUNG: Diese Anleitung muss vor dem Einsatz des Produkts aufmerksam gelesen und vollständig verstanden werden.
ATENCIÓN: Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.
ATENÇÃO: Os operadores devem ler e entender completamente este manual antes de usar o produto.
ΠΡΟΣΟΧΗ: Οι χειριστές αυτού του προϊόντος πρέπει να διαβάσουν και να καταλάβουν πλήρως τις οδηγίες του εγχειριδίου πριν από την χρήση του.
الحذرن: على العمال قراءة وفهم هذا الدليل بكامله قبل البدء باستعمال المنتج.

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FEATURES

The ear washing syringes GIMA are instruments used for ORL irrigation for ear washing and withdrawal of earwax plugging. Models of different capacities are available, completely made of metal or with glass bowl with graduated scale. All models are equipped with cannulas with interchangeable screw, in order to be better suitable to all typologies of use.

PRESCRIPTIONS

The product must be used by specialized medical personnel only, adequately trained to use and service this kind of products.

Particular attention must be paid as regards the instruments care. The incorrect use or the lack of proper maintenance can make the instruments quickly deteriorate.

The product must be sterilized before using it. It is recommended to sterilize it in autoclave.

UNPACKING

Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled.

They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

OPERATION

First of all, it is important to choose the most suitable nozzle for the required use and screw it down into the appropriate hole. If the corrected nozzle is already installed onto the instrument, check whether it has been thoroughly screwed down.

Use the instrument only for irrigation or for the intended uses!

An improper use or by unqualified personnel can cause injuries to the patient!

After using it, unscrew the nozzle, clean and sterilize the instrument.

MAINTENANCE

It is important to carefully examine the product before using it to monitor the presence of breaks, defects or wear and tear. Particular attention must be paid to moving parts, tips and connectors. The instruments must be kept in appropriate boxes or protected trays, in dry and clean environment.

In case the product is defective or excessively worn out, do not use it and do not try to repair it.

For repair apply to your own retailer, GIMA or qualified personnel.

After using it, the product must be thoroughly cleaned with warm water, disinfected and rinsed possibly with distilled water, carefully dried and sterilized before the next use.

Be careful that the product does not get in contact with acids or aggressive disinfectants, which can corrode it.

We guarantee the product for a minimum of 20 sterilizing cycles, if they are sterilized in compliance with the following criteria. Most instruments fully pass this number of cycles in presence of a care and maintenance protocol.

The recommended sterilizing method is steam sterilization in autoclave; the cycle temperature must not exceed 134°C (275°F) for up to 10 minutes, to avoid damaging the product.

It is recommended to follow the instructions of the autoclave manufacturer and arrange the instruments in such a way that their surface is directly and fully exposed to steam. Make them cool before using it. Hot air sterilization is not recommended.

GUARANTEE

Thank you for having acquired our product. This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, etc. No compensation will be given for loss of use of the product.

Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.

The guarantee will be invalidated if the registration number has been removed, cancelled or altered. Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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