# seca

## seca 354

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### Warranty

A two-year warranty from date of delivery applies to defects attributable to poor materials or workmanship. All moveable parts - batteries, cables, mains units, rechargeable batteries etc. - are excluded. Defects which come under warranty will be made good for the customer at no charge on production of the receipt. No further claims can be entertained. The costs of transport in both directions will be borne by the customer should the equipment be located anywhere other than the customer's premises. In the event of transport damage, claims under warranty can be honoured only if the complete original packaging was used for any transport and the scale secured and attached in that packaging just as it was when originally packed. All the packaging should therefore be retained.

A claim under warranty will not be honoured if the equipment is opened by persons not expressly authorised by seca to do so. We would ask our customers abroad to contact their local sales agent in the event of a warranty matter.

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seca operates worldwide with headquarters in Germany and branches in:  
**seca france**  
**seca united kingdom**  
**seca north america**  
**seca schweiz**  
**seca zhong guo**

**seca nihon**  
**seca mexico**  
**seca austria**  
**seca polska**  
**seca middle east**  
**seca suomi**  
**seca américa latina**  
**seca asia pacific**

**seca denmark**  
**seca benelux**  
**seca lietuva** and with exclusive partners in more than 110 countries.

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### Declaration of conformity

seca gmbh & co. kg hereby declares that the product meets the terms of the applicable European directives. The unabridged declaration of conformity can be found at: [www.seca.com](http://www.seca.com).

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




## Medical Measuring Systems and Scales since 1840

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## Set up scale

- Please remove all packaging and insert the batteries as described below.
- Place the scale on a level and stable surface. Then position the tray as shown on the base of the scale and pull forward until the tray is heard to snap into place, and the release button is in the top position (Fig. 1). The tray is now securely interlocked with the base.
- To remove the tray just press the release button hard and slide back the tray (Fig. 2).

## Connect power supply

Batteries are used to supply power to the scale. The 4 Mignon (AA) batteries supplied should last for up to 20,000 weighing operations. When inserting batteries, always ensure the correct polarity!

- Take off the tray and then turn the scale over carefully so that the base is accessible.
- Open the battery compartment lid and insert the batteries supplied according to the polarity shown on the battery compartment (Fig. 3).
- Now close the battery compartment cover again carefully.

## Correct weighing procedure

The scale is equipped with an "AUTOHOLD" function to ensure correct readings for weights exceeding 500g. This function provides for automatic "fixing" of the weighing result on the digital display until the scale switches off after 2 minutes or is used for another weighing operation.

- Switch on the scale by pressing the green key. "0" should now appear in display. The signal tone indicates that the scale is now ready for weighing.
- When weighing infants: Place the infant on the tray, ensuring that the child lies as still as possible (Fig. 4).
- When weighing children able to stand up: Remove the tray as described. Check the secure positioning of the base. Place the child on the scale and make sure he/she keeps still (Fig. 5).
- Read the weighing result off display. The scale must not be subjected to weights exceeding the maximum load specified!
- The integrated TARE function can be used to tare in items present on the scale such as bedding, nappies etc. Their weight will then be automatically deducted from the total weight measured:
  - First place the item to be tared in on the scale.
  - Then press key **TARE**. The display of the scale is now reset to "0" and a normal weighing operation can begin.
- **Note:** The accuracy of the differential measurement may be impaired if the infant is kicking vigorously.
- The scale can be switched off by pressing the **START** button.

## Use of differential display function

- After finishing the weighing process, press the **STO** until the indicator shows -> **M**, and a signal tone is heard. The weight value is now saved and will also remain stored once the scale is switched off.
- With a subsequent weighing operation it is then possible to switch the display over to differential display, i.e. display of current value less saved value, by briefly pressing the **STO** button. You can use this function for example to measure the intake of milk when breastfeeding or have Baby's increase in weight displayed since being last weighed.

## Automatic switch-off/energy-saving mode

The scale remains ready for weighing for approx. 20 seconds after being switched on. Once weighing has been performed, the value is displayed for approx. 2 minutes until the scale switches off automatically.

## Switchover of weight display

It is possible to switch over the weight display between kilograms (kg) and pounds: ounces (lbs:oz) and decimal pounds (lbs) using the **kg / lbs** button.

## Cleaning

Clean the tray and the outside of the scale with a damp cloth as required. Under no circumstances use abrasive or acid cleaning materials, white spirit, benzene or similar substances for cleaning!

## What to do if ...

- ... **no weight is displayed under load?**
  - Check that the batteries are inserted correctly.
- ... **"batt" is displayed?**
  - The battery voltage will not last much longer. Insert a new set of batteries.
- ... **the display "Stop" appears?**
  - Maximum load has been exceeded!
  - Remove the weight on the scale as quickly as possible in order to avoid permanent damage to the scale.

## Disposing of batteries

Dispose of batteries via local collection points, not with domestic waste. When purchasing new batteries, select those low in harmful substances which do not contain mercury (Hg), cadmium (Cd) or lead (Pb).

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